



THE WELSH PONY & COB SOCIETY OF NEW ZEALAND (INC)
DNA COLLECTION SHEET

FOR LAB USE ONLY

Case #

Date Received:

Name of Animal

Reg No

Sex

Colour

Date of Birth

Name of Sire

Reg No

Name of Dam

Reg No

Fold Here

SAMPLE SECTION

PLACE HAIR
ROOTS HERE
(20-30 HAIRS)

TAPE
CENTRE OF
HAIR
SHAFTS
HERE

Fold Here

Instruction for the collection of hair sample for DNA TYPING

1. One envelope for each horse
2. Pull mane hair from above the withers (about mid neck region).
3. Firmly grasp a tuft of 25/40 hairs as close to the skin as possible
4. Pull evenly and directly up and away from the neck so that the hair will tend to come out by the roots rather than break off.
5. Look closely at the ends to ensure the hair roots are resent, use a magnifying glass if necessary. At least half the hair sample should have roots
6. **Avoid contaminating with hair from other horses**
7. Collected hair should be **immediately placed** in labeled envelope and sealed (avoid licking)

